



Email : sales@laboratoryinstrumentindia.com

Product Name :
Centrifugal Pumps Study Unit

Product Code :
LBNY-0005-16000012



www.laboratoryinstrumentindia.com

Description :

Equipped with all the instruments and controls necessary for testing, it requires only to be connected to electrical supply and filled with clean water. The trolley-mounted self-sufficient unit has been designed to study the features of centrifugal pumps.

The group is supplied with manuals which describe all parts of the unit, the installation and utilisation procedures, as well as many exercises with the relative results.

The Centrifugal Pumps Study Unit is supplied in four different versions, the first two of which allow also for series-parallel tests:

One pump with a.c. motor

One pump with d.c. motor

Two pumps with 2-speed a.c. motors

Two pumps with one 2-speed a.c. motor and one d.c. motor

Technical Specification :

Centrifugal Pumps Study Unit



Laboratory Instrument India

