



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Diaphragm And Piston Pumps Test Unit

**Product Code :**  
LBNY-0005-16000015



**Description :**

The Diaphragm And Piston Pumps Test Unit is completely self-sufficient and it is fitted with a wheeled frame the working liquid is the water. The reciprocating pumps fitted with the bench are micro metrically regulated by the piston upstroke this allows to change the volume of the pumping chamber and consequently to modify the delivered capacity without reducing the delivery pressure.

The group is supplied with manuals which describe all parts of the Diaphragm And Piston Pumps Test Unit, the installation and utilisation procedures, aswell as many exercises with the relative results. The Diaphragm And Piston Pumps Study unit allows the student to understand and to experiment the reciprocating pumps running.

**Technical Specification :**

Weight and dimensions:

Dimensions: 550 x 750 x 1500 h mm

Weight: 130 kg.



**Laboratory Instrument India**