



Email : sales@laboratoryinstrumentindia.com

Product Name :
Gear Pumps Study Unit

Product Code :
LBNY-0005-16000016



Description :

The group is supplied with manuals which describe all parts of the Gear Pumps Study Unit the installation and utilisation procedures, as well as many exercises with relative results. The gear pump mounted onto the bench is characterized by the possibility of regulating the delivered flow rate, keeping constant the delivery pressure, changing the rotation speed of the operation electrical motor. The gear pump study unit allows the student to understand and to experiment the oil gear pumps running.

Technical Specification :

Required services:

Electrical supply: 380 V AC, three-phase, 50/60 Hz, 1.2 kW.

Weight and dimensions:

Dimensions: 1100 x 550 x 1750 h mm

Weight: 110 kg.



Laboratory Instrument India