



Email : sales@laboratoryinstrumentindia.com

Product Name :
Interface-Portable

Product Code :
LBNY-0004-780074



Description :

Interface-Portable

Technical Specification :

Interface-Portable
Full color touchscreen: 10 inch
Frequency: 1.2 GHz quad-core processor
Memory: 1.5 GB RAM, 16 GB
Built-in Wi-Fi
Built-in Bluetooth,
USB connections
Can export data to third-party apps
Suitable for wet and dry labs
Operates online or offline
Software is included
Lab Manager classroom management application included
Front & back cameras
Built-in GPS, Accelerometer, Microphone, & Speaker
Port for 2 digital sensors
Can simultaneously connects to 8 or more wireless sensors.



Laboratory Instrument India