



Email : sales@laboratoryinstrumentindia.com

Product Name :
Head Loss Study Unit

Product Code :
LBNY-0005-16000018



www.laboratoryinstrumentindia.com

Description :

The Head Loss Study Unit is fitted with contraction devices for the determination of flow-rate, a Pitot tube to measure the velocity of the fluid flow, and one manometers to measure pressure differences. The Head Loss Study Unit combines all the elements necessary to conduct extensive studies on head losses in pipes and valves. It is easy to handle and has different test shapes and valves of different kinds. The measures can be made in different operating conditions and the values obtained can be easily compared with theoretical results.

The training unit is supplied in two versions:

Fed from water network

Self-sufficient

Technical Specification :

Required service:

Water supply: 3000 lt/h.

Water discharge: rubber holder \hat{A} $\frac{3}{4}$ min.

Electrical power supply: 220 V a.c. single phase, 50/60 Hz

Weight and dimensions:

Dimensions: 2650 x 900 x 1900 h mm

Net weight: 140 kg.



Laboratory Instrument India