



Email : sales@laboratoryinstrumentindia.com

Product Name :
Apparatus For Flow Rate Measurements

Product Code :
LBNY-0005-16000019



www.laboratoryinstrumentindia.com

Description :

The apparatus consists of a hydraulic feed bench on which some instruments in series are tested and it allows the measuring of the same flow rate verifying in this way, the imprecision degree of flow rate measuring methods. The Apparatus For Flow Rate Measurements has been designed to allow the comparison between different methods of flow rate measurements. The group is supplied with manuals which describe all parts of the Apparatus For Flow Rate Measurements unit, the installation and utilisation procedures, as well as many exercises with the relative results. It is a useful educational instrument because it allows the student to acquire the fundamental knowledge concerning the hydraulic measurements; the student can also realize by himself about the uncertainty of the measurement.

Technical Specification :

Service requirements:

Water supply: before the beginning of the test.

Electrical power supply: 220/240 V single-phase, 50/60 Hz, 0.5 kW.

Weight and dimensions:

Hydraulic feed bench dimensions: 2000x710x1100 h mm

Net weight: 100 Kg.

Laboratory Instrument India

