



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Flow Visualization Mini Channel

**Product Code :**  
LBNY-0005-16000021



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

It can be fed from our Hydraulic Bench or from the mains and it allows, by making use of different weirs, to carry out experiments on the open channels. The Flow Visualization Mini Channel unit, made of transparent plexiglas, allows to study the behaviour of a water flow in a channel of small dimensions. An ink jet device allows to visualize the water flow around the models.

**Technical Specification :**

**Services required:**

Water supply from mains: 300 l/h.

Weight and dimensions:

Dimensions: 630x15x400 h mm

Weight: 15 kg.



**Laboratory Instrument India**