



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Apparatus For The Determination Of Pressure Losses In Pipes

**Product Code :**  
LBNY-0005-16000022



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Hydraulic feed is obtained through an adjustable height reservoir at constant hydraulic head. Pressure losses are measured by means of two U manometers, one water type (800 mm) and the other mercury type (400mm). The Apparatus For The Determination Of Pressure Losses In Pipes has been designed to determine the pressure loss in a pipe in which water is made to flow.

**Experiments:**

Measurement of pressure loss in a pipe as a function of flow-rate.  
Change from lamnar to turbulent flow.  
Verification of Reynolds number.

**Technical Specification :**

**Required services:**

Hydraulic supply: 100 l/h from mains.  
Weight and Dimensions:  
Dimensions: 700 x 400 x 1700 h mm.  
Weight: 15 kg.

---

# Laboratory Instrument India

