



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Bluetooth Adapter

**Product Code :**  
LBNY-0004-780075



**Description :**

Bluetooth Adapter

**Technical Specification :**

Bluetooth Adapter:

Adapter is compatible with the wireless sensors in order to have connection from sensors to PC or laptop.

Start Using Bluetooth—Adds Bluetooth to your Windows 11, 10, 8, and 7 PC that doesn't already have it.

Connect Your Devices—Supports wireless Bluetooth headphones, speakers, headsets, keyboards, mice, phones, and more.

Tiny Design—This mini USB dongle can stay plugged in while you are on the move. This adapter is intended only for computers that do not already have built-in Bluetooth.

Compatibility—Not usable in car stereos, TVs, or Mac computers. Not recommended for employer-provided computers. Does not add Bluetooth 4.0 support to Windows 7. Special drivers (provided) are needed in Windows 7 for Bose and some other recent devices.



**Laboratory Instrument India**

