



Email : sales@laboratoryinstrumentindia.com

Product Name :
Spirometer Digital Sensor

Product Code :
LBNY-0004-780076



Description :

Spirometer Digital Sensor

Technical Specification :

Spirometer Digital sensor
Measures airflow, duration and lung volume
Supplied with disposable mouth pieces
Sample rate range between 50 to 100 Hz
Read the volume in liters
Bi-directional airflow
Supplied with a pack of 100 mouth pieces.

Features:

Displays volume in liters
Minimal resistance to airflow
Bi-directional airflow (inspiration and expiration)

Applications:

Conduct respiratory experiments
Compare a student's airflow before and after exercise
Investigate the lung volume of athletes vs. non-athletes

Compare musicians and non-musicians
Determine total lung capacity.



Laboratory Instrument India