



Email : sales@laboratoryinstrumentindia.com

Product Name :
Oxygen O2 Gas Digital Sensor

Product Code :
LBNY-0004-780077



Description :

Oxygen O2 Gas Digital Sensor

Technical Specification :

Oxygen Gas Digital Sensor

The O2 Gas Sensor measures oxygen concentration in air, and is temperature compensated for use in different environments.

It has a wide measurement range, allowing it to be used to study human and cellular respiration.

The O2 Gas Sensor also includes a 250 mL bottle that can be used to study the rusting of iron or as a respiration chamber for monitoring plants and insects.

Measures gas concentration in open or closed system

Range between 0 to 100%

Accuracy: +/- 1%

Resolution: 0.025%

Supplied with extension cable.



Laboratory Instrument India

