



Email : sales@laboratoryinstrumentindia.com

Product Name :
Yaw Probe

Product Code :
LBNY-0005-1620006



www.laboratoryinstrumentindia.com

Description :

Featured machined aluminum triangle that shows the orientation of the ports from outside the test Section and acts as a handle

Three ports

Mounts easily in the test section of an Educational Wind Tunnel

Length of 10.5 inches (26.7cm)

Outside diameter (O.D.) of 0.125 inches (3.175mm)

Comes with ample tubing.

Technical Specification :

Yaw Probe



Laboratory Instrument India