



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Yaw Probe

**Product Code :**  
LBNY-0005-1620006



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Featured machined aluminum triangle that shows the orientation of the ports from outside the test Section and acts as a handle

Three ports

Mounts easily in the test section of an Educational Wind Tunnel

Length of 10.5 inches (26.7cm)

Outside diameter (O.D.) of 0.125 inches (3.175mm)

Comes with ample tubing.

**Technical Specification :**

Yaw Probe



**Laboratory Instrument India**