



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pressure Wing

**Product Code :**  
LBNY-0005-1620008



**Description :**

18 Flush-mounted pressure taps around its Clark Y-14 airfoil  
Angular position is adjustable  
Mounts vertically and spans the test section of an Educational Wind Tunnel (12 inches / 30.48cm)  
Comes with ample tubing  
Pressure Wing has a chord of 3.5 inches (8.9cm).

**Technical Specification :**

Pressure Wing



**Laboratory Instrument India**