



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pitot Static Probe

**Product Code :**  
LBNY-0005-1620009



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Features a rounded-tip total-pressure tap  
Reaches 2.6 inches (6.6cm) forward of the bend  
Has an OD of 0.125 inches (3.175mm)  
6 static ports  
Extends 13 inches (33cm) to the bend  
Comes with ample tubing.

**Technical Specification :**

Pitot Static Probe



**Laboratory Instrument India**