



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pitot Static Probe

Product Code :
LBNY-0005-1620009



www.laboratoryinstrumentindia.com

Description :

Features a rounded-tip total-pressure tap
Reaches 2.6 inches (6.6cm) forward of the bend
Has an OD of 0.125 inches (3.175mm)
6 static ports
Extends 13 inches (33cm) to the bend
Comes with ample tubing.

Technical Specification :

Pitot Static Probe



Laboratory Instrument India