



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wireless Light Sensor

Product Code :
LBNY-0004-780078



Description :

Wireless Light Sensor

Technical Specification :

Wireless light sensor

The Wireless Light Sensor features two separate apertures – one for ambient light measurements and one for directional light measurements.

The ambient sensor measures illuminance and UV Index, while the spot (directional) aperture measures light level and color intensity. Our software displays the relative intensities of Red, Green, and Blue light, then sums them to determine the level of White light.

Photodiode spectral response from 300 to 1100nm for normal light

Can also measure: RGB, UVA, UVB and UV index

Range: 0 to around 130,000 LUX



Laboratory Instrument India