



Email : sales@laboratoryinstrumentindia.com

Product Name :
Transonic Wind Tunnels

Product Code :
LBNY-0005-16200018



Description :

Employing variable test section porosity, computer-controlled valves, chokes and flaps, target airspeeds are easy to achieve and maintain. Using compressed air Transonic Wind Tunnels offer running speeds of $M=0.3$ to $M=1.3$. Systems with test section sizes of 8-inch x 8-inch (20.32cm x 20.32cm) and 20-inch x 20-inch (51cm x 51cm) are available.

Technical Specification :

Nominal Mach Number Range: $M=0.3$ to $M=1.3$.



Laboratory Instrument India