



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Set of Electricity Sensors

**Product Code :**  
LBNY-0004-780079



**Description :**

Set of Electricity Sensors

**Technical Specification :**

Set of Electricity sensors  
Wireless Current sensor  
Two ranges: +/- 1 A and +/- 0.1 A or better  
Resolution 0.2 mA (1 A) and 0.02 ( 0.1A)  
Includes 1 red and 1 black shrouded banana-to-alligator test leads  
Wireless Voltage sensor  
Two ranges: +/- 15 V and +/- 5V  
Resolution: 7mV (15 V) and 2mV  
Includes 1 red and 1 black shrouded banana-to-alligator test leads.



**Laboratory Instrument India**