



Email : sales@laboratoryinstrumentindia.com

Product Name :
Set of Electricity Sensors

Product Code :
LBNY-0004-780079



Description :

Set of Electricity Sensors

Technical Specification :

Set of Electricity sensors
Wireless Current sensor
Two ranges: +/- 1 A and +/- 0.1 A or better
Resolution 0.2 mA (1 A) and 0.02 (0.1A)
Includes 1 red and 1 black shrouded banana-to-alligator test leads
Wireless Voltage sensor
Two ranges: +/- 15 V and +/- 5V
Resolution: 7mV (15 V) and 2mV
Includes 1 red and 1 black shrouded banana-to-alligator test leads.



Laboratory Instrument India