



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Wireless Temperature Sensor

**Product Code :**  
LBNY-0004-780080



**Description :**

Wireless Temperature Sensor

**Technical Specification :**

Wireless Temperature Sensor  
Stainless steel temperature sensor  
Range : -40 °C to +125 °C, at least  
Accuracy: +/- 0.5 °C  
Resolution : 0.01 °C

**Features:**

Simply pair and go, no cables or adapters to manage.  
Variable sampling rate for capturing small, fast changes or experiments that run for hours, days, or weeks.  
Bluetooth connectivity and long-lasting coin cell battery.  
Logs temperature data directly onto the sensor for long-term experiments.  
Dust, dirt, and sand-proof and water resistant.



**Laboratory Instrument India**

