



Email : sales@laboratoryinstrumentindia.com

Product Name :
CCTV Tester

Product Code :
LBNY-0002-610011



Description :

CCTV Tester

Technical Specification :

Quickly measure the length of power, BNC cable, network cable, telephone cable, RVV cable and elevator cable, cat 5/6 cable's over a distance of up to 1200 meters.

Detects cable whether it is unblocked and locates breakpoints precisely.

With the attached cable tester tool, you can test the arrangement type of the network cable and whether it is unblocked.

CCTV Camera Tester, 4K IP Camera Tester,
IPS Touch Screen Support POE Wifi 4IC H.265 HDMI
8GBTF Card DC 5V/12V/24V Output IPC Tester

CCTV Tester, CCTV Tester Equipments, CCTV Tester Tools, CCTV Tester Tool Kits, CCTV Tester Manufacturers, CCTV Tester Suppliers from India, China, Kenya



Laboratory Instrument India

