



Email : sales@laboratoryinstrumentindia.com

Product Name :
Blowby Meters

Product Code :
LBNY-0005-16500013



www.laboratoryinstrumentindia.com

Description :

Blowby Meter is an ideal solution to the problem of easily measuring the gas flow out of engine crankcases as it can be incorporated into engine ventilation systems in a low intrusive manner. Engine deterioration may be assessed in several ways, one of which is to measure cylinder gas loss; the accumulation of crankcase gas is a result of this loss, and it provides an indication of wear in the valve stems, valve guides, piston rings, and cylinder bores.

Features:

- Continuous Flow Measurements
- Vortex Shedding Technology
- Very Low Pressure Drop (Flow Range 4~5650 l/min)
- Small Physical Size
- 3 1/2 Digit LCD Display
- Excellent at Low Flows
- Linearising Electronics
- High Flow Alarm
- Analogue Output (0-2V or 0-10V, selectable)
- Pulse Train Output
- 40 : 1 turndown ratio.

Technical Specification :

Blowby Meters



Laboratory Instrument India