



Email : sales@laboratoryinstrumentindia.com

Product Name :
Three Closed Cup Flash Point Apparatus

Product Code :
LBNY-0005-1660004



Description :

An automatic flash detector reduces the chance of misinterpretation of the test result by the operator, and enhances repeatability. Gas for the test jet/shutter assembly is supplied from an integral gas tank, which is filled using a standard butane (lighter) refill cartridge. The digital display shows test parameters, instrument status, and the test result. An audible prompt is sounded when an action is required from the operator. The Series 3 flash point tester is lightweight and suitable for mobile use when operated from the optional 12 volt power supply adaptor. All functions on the Series 3 flash point tester are accessed via a two button keypad and the Seta Multifunctional And Rotational Test (SMART) control feature.

Features:

- Audible and on-screen prompts
- Automatic flash detection
- Small sample size - 2 or 4ml
- Simple calibration
- Portable operation.
- Ambient to 300°C
- Flash point test in less than 2 minutes.

Technical Specification :

Three Closed Cup Flash Point Apparatus



Laboratory Instrument India