



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital Baby Scale

Product Code :
LBNY-0004-790002



Description :

Digital Baby Scale

Technical Specification :

Weighing Capacity 300g.

Easy to Read LED display.

Zero calibration.

Electronic damping facility for eliminating the reading fluctuations caused by moving baby.

Weighing pan suitable for weighing new born babies and the construction not allow the baby to slip from the tray.

Automatic switch off.

Hold the measured value irrespective of the baby movements.

Option: breast milk intake function.

Tare, Auto-hold, and kg/ibs functions.

Hospital use model.

Provide samples for the experiment.

The pan facility to measure the length of baby.

Operate on input supply 200 to 230 Vac, 50 Hz.

Laboratory Instrument India

