



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mini Turbine Unit

Product Code :
LBNY-0005-1720008



Description :

The unit consists of a pump and a storage tank on a steel base input power sensors, and an interface box. The box includes power supply to the motor and computer interface unit. Mini Turbine Service Unit service module provides a water supply and input power measurement for the mini turbine demonstration unit. It also provides, interface ports to the sensors and computer. Both are separately supplied. The Mini Turbine Service Unit is to be used with Mini turbine demonstration unit, and a computer.

Technical Specification :

Sensors with display and computer interface unit: Pressure and flow rate.
Maximum flow rate: 80 lpm.
Power supply: 220 V 1 Ph 50 Hz.



Laboratory Instrument India