




Email : sales@laboratoryinstrumentindia.com

Product Name :
Educational Cavitation Tunnel Instrument

Product Code :
LBNY-0005-1790004



Description :

High Speed / Frame Rate Camera system for observing cavitation.
Window mounting hydrophone system for use in the educational cavitation tunnels.
Educational Wake Rake with integrated pressure transducers for measuring the wake profile within the tunnel.
Flat Plate and Cylinder Model Experiment for use in Tunnels.
Pitot Static tube for use in measuring section of Tunnel.

Technical Specification :

Educational Cavitation Tunnel Instrument



Laboratory Instrument India