



Email : sales@laboratoryinstrumentindia.com

Product Name :
Level Process Training System

Product Code :
LBNY-0005-1850006



Description :

The Level Process Training System is a compact and mobile unit for a wide range of experiments in flow control. It gives students a greater understanding of the stability of simple control systems. The self-contained unit can do many experiments, but it can also connect to other products in the range for extra experiments. For cascade control of flow and level, it can link to the optional Flow Process Training System. For distributed control, it can connect to the optional Computer Control System

Technical Specification :

The main parts are:

- Industrial controller with auto-tune feature
- Current-to-pressure converter
- Differential pressure transmitter
- Two-channel chart recorder
- Pneumatic control valve
- Three-speed pump
- Transparent vessel
- Reservoir.

Laboratory Instrument India

