



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital Interface Trainer

Product Code :
LBNY-0005-18500014



Description :

This training kit contains pulse generator of different frequency levels, DC power supply system, data switches and breadboard. Ideal for projects of schools and colleges, this trainer has inbuilt PC board for its trouble free operation. Offered Digital Interface Trainer acts as unique computer system. It serves as a wonderful medium to learn how to write into ROMs and RAMS. Its digital section contains buffered LED readout system, logic switches and clock frequency system. We are an eminent supplier and manufacturer of best grade Digital Interface Trainer. The entire structure of this trainer is installed on metal base that has locks and handles.

Features:

- Error free operation
- Single PC board based design
- Fast functioning.

Technical Specification :

Digital Interface Trainer



Laboratory Instrument India

