




Email : sales@laboratoryinstrumentindia.com

Product Name :
Engine Speed Control Apparatus

Product Code :
LBNY-0005-18500020



Description :

Experiments:

The measurement of system dynamics from step response information.

Inner loop feedback compensation.

The use of dither signals in the compensation of system non-linearities.

P+I controller design.

It is good for extended or advanced control experiments, and is ideal for student project work. The flexible design of the equipment allows the user to develop many other analysis and control exercises to suit their needs.

Technical Specification :

Engine Speed Control Apparatus



Laboratory Instrument India