



Email : sales@laboratoryinstrumentindia.com

Product Name :
Coupled Drives Apparatus

Product Code :
LBNY-0005-18500021



Description :

Experiments:

Simultaneous control of speed and tension

Independent control of speed and tension

The flexible design of the equipment allows the user to develop many other analysis and control exercises to suit their needs. It is good for extended or advanced control experiments, and is ideal for student project work.

Practical methods of controlling multi-variable electro- mechanical systems.

Technical Specification :

Coupled Drives Apparatus



Laboratory Instrument India