



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mobile Hydraulics

Product Code :
LBNY-0005-1860003



Description :

The program includes: Instructional Modules, Audio Nisual Presentations, General Reference Manuals, Instructor's Guide, Hydraulic Tester and Mobile Hydraulic Workstation. The Mobile Hydraulics provides hands on activities with components commonly used in construction equipment, agricultural and utility vehicles. This self contained unit is a turnkey program with all the various components, hoses, pumps, and trays to build dozens of circuits. It has a variable working pressure capacity of 500-1000 and includes a heavy duty workbench with instruments and gauges.

Technical Specification :

Mobile Hydraulics



Laboratory Instrument India