



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Truck Anti lock Brakes

**Product Code :**  
LBNY-0005-1860004



**Description :**

Truck Antilock Brakes is the latest Mentor setup fully functional with all sensors and fault insertion. It follows the factory training program. A power supply and fault panel is located on the rear of the unit. Truck Anti lock Brakes the display is the shadow box type, with all air brake components on the front and power supply, with storage cabinet located on the rear. The rear of the display is painted with melamine paint. Silk screens and actual reluctors convey all the technology in this working trainer. The cabinet main structure is cabinet grade plywood covered with Formica on the front, sides, and top.

**Technical Specification :**

Truck Anti lock Brakes



**Laboratory Instrument India**