



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Engine Test Stand With CNG

**Product Code :**  
LBNY-0005-18600010



**Description :**

The stand is to accommodate radiator mount with varied positioning slots and sliding brackets positioned in the front of the engine test stand. The engine mount is a heavy duty, test stand, constructed from 3" x 3" steel mechanical tubing frame, 11 gauge heavy wall type, welded at all joints. The entire test stand is positioned on specially designed vibration dampers and are attached to the main frame.

**Technical Specification :**

Engine Test Stand With CNG



**Laboratory Instrument India**