



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Automotive Solar Ventilation Trainer

**Product Code :**  
LBNY-0005-18600015



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The module includes the components as used in the Prius scaled down for quick and easy instruction. Included is a 12-20 Volt Solar Panel, Solar ECU, Blower, Relay High Output Blower Fan and Ventilation Control Switch. Automotive Solar Ventilation Trainer /program is designed to demonstrate the basic operation of the optional Solar HVAC assist offered. The air inlet/outlet Servo function is simulated with colored lights to show timed operation while the actual air flow comes through the live dash vent.

**Technical Specification :**

Automotive Solar Ventilation Trainer



**Laboratory Instrument India**