



Email : sales@laboratoryinstrumentindia.com

Product Name :
Automotive Solar Ventilation Trainer

Product Code :
LBNY-0005-18600015



www.laboratoryinstrumentindia.com

Description :

The module includes the components as used in the Prius scaled down for quick and easy instruction. Included is a 12-20 Volt Solar Panel, Solar ECU, Blower, Relay High Output Blower Fan and Ventilation Control Switch. Automotive Solar Ventilation Trainer /program is designed to demonstrate the basic operation of the optional Solar HVAC assist offered. The air inlet/outlet Servo function is simulated with colored lights to show timed operation while the actual air flow comes through the live dash vent.

Technical Specification :

Automotive Solar Ventilation Trainer



Laboratory Instrument India