



Email : sales@laboratoryinstrumentindia.com

Product Name :
Manual Air Condition

Product Code :
LBNY-0005-18600020



www.laboratoryinstrumentindia.com

Description :

The components of this unit are pre-assembled with a compressor operating with refrigerant. A large electric motor is mounted with a control housing in a self contained heavy duty roller cart. Manual Air Conditioning Trainer allows the student to learn the basic test and maintenance procedure for the refrigerant loop. Charging checking for leaks, A replacing components, servicing parts, and other tasks are covered in this program. Included are a set of charging and manifold gauges and battery cables.

Technical Specification :

Manual Air Condition



Laboratory Instrument India