



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Manual Air Condition

**Product Code :**  
LBNY-0005-18600020



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The components of this unit are pre-assembled with a compressor operating with refrigerant. A large electric motor is mounted with a control housing in a self contained heavy duty roller cart. Manual Air Conditioning Trainer allows the student to learn the basic test and maintenance procedure for the refrigerant loop. Charging checking for leaks, A replacing components, servicing parts, and other tasks are covered in this program. Included are a set of charging and manifold gauges and battery cables.

**Technical Specification :**

Manual Air Condition



**Laboratory Instrument India**