



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electronic Climate Control AC Trainer

Product Code :
LBNY-0005-18600022



Description :

Electronic Climate Control AC Trainer is an automotive air conditioning program designed to offer hands on air conditioning systems for students. The trainer consist of an electronic climate control air conditioning unit including: Compressor, blower, fan, evaporator condenser, accumulator, pressure switch and related parts mounted on an industrial duty mobile cart. The components are based on an AC unit, complete with electronic controls, and comes precharged with refrigerant. Included are a set of charging hoses and manifold gauges and battery cable.

Technical Specification :

Electronic Climate Control AC Trainer



Laboratory Instrument India