



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Automatic Transaxle Trainer

**Product Code :**  
LBNY-0005-1860023



**Description :**

The trainer features a fully functioning Automatic Transaxle Trainer that can be manually shifted through all foreword and reverse gears. It features the popular Automatic Transaxle found on most chrysler front wheel drive cars. The Automatic Transaxle Trainer is designed to assist in educating the fundamentals of a front drive transaxle. The package includes a complete competency based training program which instructs the student from basic operation and adjustment to complete overhaul of the transmission.

**Technical Specification :**

Automatic Transaxle Trainer



**Laboratory Instrument India**