




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hybrid Vehicle Demonstrator

**Product Code :**  
LBNY-0005-18600024



**Description :**

The Hybrid Vehicle Demonstrator includes a Hybrid engine, a battery and an Electronic Control Unit, several meters showing power flow between the units, a brake switch and a potentiometer mimicking the accelerator pedal. The system can also be used to show the effects of regenerative braking. Two models are available: a bench top version and a panel mounted version. The controls the system to show students how the power is routed in a hybrid depending on the State of Charge of the Battery, and to allow them to make measurements on the engine performance under different load conditions.

**Technical Specification :**

Hybrid Vehicle Demonstrator



**Laboratory Instrument India**