



Email : sales@laboratoryinstrumentindia.com

Product Name :
Engine Repair Program

Product Code :
LBNY-0005-18600025



www.laboratoryinstrumentindia.com

Description :

The engine stand provides a stable work bench for major engine service suitable for use by schools, professional mechanics and engine re-builder. Engine Repair Program provides hands on training for those who wish to specialize in engine repair. The engine stand features a rugged tubular frame which supports up to 1,250lbs Include engine, stand, and tools for overhaul with curriculum. The stand has 250 degree rotation heads which lock in any position, a 23" handle to allow easy rotation to the engine and 3" diameter steel wheels to provide smooth maneuvering.

Technical Specification :

Engine Repair Program



Laboratory Instrument India