



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Automotive Electro Trainer

**Product Code :**  
LBNY-0005-18600045



**Description :**

Automotive Electro Trainer provides a "hands-on" working knowledge of voltage, current, and resistance with series and parallel circuits. Over 250 of these currently used to train. Automotive Electro Trainer kit allows students to gain a foundation in the principles of electricity and electronics. Graphic symbols are screened into each component's acrylic background. This helps students to understand circuit reading and wiring diagrams. It provides a "hands on" working knowledge of voltage, current, and resistance with series and parallel circuits.

**Technical Specification :**

Automotive Electro Trainer



**Laboratory Instrument India**