




Email : sales@laboratoryinstrumentindia.com

Product Name :
Classic Series Trainers

Product Code :
LBNY-0005-18600048



Description :

The courseware is competency based including; information sheets, instruction sheets, student self-check, instructor's final checklist, and check out activities. Classic Series Trainers is made from live components, mounted on a Mylar (erasable) surface, with heavy rubber trim and are predrilled for steel wall mount brackets for classroom security. Each board can keep 2-3 students active doing hands on skill development.

Technical Specification :

Classic Series Trainers



Laboratory Instrument India