



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Air Bag Trainer Redeployable

**Product Code :**  
LBNY-0005-18600049



**Description :**

With the Air Bag Trainer, the student can learn the operation, testing and troubleshooting of an actual operating deployable air bag.

**Features:**

Crash and safety sensors are directly available on the board, without the student getting involved with a live vehicle. The teachings of "Back-Up" power supply are shown on the board with relative ease, in case of battery destruction. Pneumatic operated deployable air bag has a flow rate valve that controls the speed of deployment into a safe range of operation. Components are displayed and accessible. The unit provides for an Air Bag indicator light that flashes trouble codes. Diagnostic Monitoring is performed with a microprocessor which checks the wiring and electrical components.

**Technical Specification :**

Air Bag Trainer Redeployable



**Laboratory Instrument India**

