




Email : sales@laboratoryinstrumentindia.com

Product Name :
CAN Bus System And Operation Trainer

Product Code :
LBNY-0005-18600052



Description :

The solution includes component carriers, base board, a power supply, and storage trays and is supplied with a full curriculum pack. CAN Bus System And Operation Trainer kit allows a fully functioning CAN Bus System, mimicking vehicle operation, to be set up using 4 Electronics Control Units representing Instrument panel, Front ECU, Power train control, and Rear ECU. Students are tasked with setting up a fully working CAN Bus System, inserting faults and using hardware and software tools to understand fault diagnosis procedures and practice. A fifth is used for system diagnosis, releasing faults and viewing CAN bus messages.

Technical Specification :

CAN Bus System And Operation Trainer



Laboratory Instrument India