



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hybrid Demonstrator Trainer

Product Code :
LBNY-0005-18600054



Description :

The Hybrid Demonstrator Trainer includes a Hybrid engine, a battery and an Electronic Control Unit, several meters showing power flow between the units, a brake switch and a potentiometer mimicking the accelerator pedal. The system can also be used to show the effects of regenerative braking. Two models are available: a bench top version and a panel mounted version. The controls the system to show students how the power is routed in a hybrid depending on the State of Charge of the Battery, and to allow them to make measurements on the engine performance under different load conditions.

Technical Specification :

Hybrid Demonstrator Trainer



Laboratory Instrument India