



Email : sales@laboratoryinstrumentindia.com

Product Name :
Recycled Gas Engine Performance Program

Product Code :
LBNY-0005-18600058



7

Description :

Programs are designed to give practical troubleshooting experience with domestic and foreign car manufacturer's engines. With the Variable Fault Feature, the instructor can use the switches to demonstrate the loss of circuits while the engine is running. Today's engines have sophisticated microprocessor based algorithms that control engine performance. Because of these technological changes, the automotive technician must understand each subsystem's unique role in operation as well as its interdependence with other components.

Technical Specification :

Recycled Gas Engine Performance Program



Laboratory Instrument India