



Email : sales@laboratoryinstrumentindia.com

Product Name :
Fuel Injection Trainer

Product Code :
LBNY-0005-18600059



Description :

The Four-Cylinder Engine trainer should consist of a complete, fully operational Multi Point EFI system. A color coded schematic chart, under a lexon polycarbonate protective screen, of the Fuel Injection Trainer system and circuit's, should also be mounted onto the trainer for reference. Original components should be utilized and clearly laid out on a display board that is mounted on an integrated floor stand, made from powder coated electro galvanized steel, mounted onto wheels for mobility.

Fuel Injection Trainer engine sensors have a signal generator to simulate various engine conditions and a corresponding electrical meter that continuously monitors the signal. The Electrical wiring diagram is provided in a special design panel to incorporate Input / Output connections from sensors and ECU. The effect of altering the sensor condition should be immediately apparent on the meters. Open source and programmable and have the ability to connect to a PC, for programming and displaying system data.

Technical Specification :

Fuel Injection Trainer



Laboratory Instrument India

