



Email : sales@laboratoryinstrumentindia.com

Product Name :
Ref AC Control And Fault Trainer

Product Code :
LBNY-0005-1880004



www.laboratoryinstrumentindia.com

Description :

Designed for single phase operation the demonstration case is divided into two main boards:-
Control board with programmable controller and normal Run/Fault switches and indicator lamps.
Mimic Diagram with Cooling and Heating sections and room temperature and humidity potentiometers.

Technical Specification :

The following main functions and components are simulated:-

- Crankcase Heater
- High Pressure Switch
- Reversing Valve
- Refrigeration/Cooling
- Heating/Heat Pump
- Compressor
- Condensing Unit
- Expansion Device (outdoor use)
- Non return Valve (outdoor use)
- Reversible Filter drier
- Sight Glass (Heating)
- Expansion Device (indoor use)
- Non Return Valve (indoor use)
- Suction Line Accumulator
- Defrost Facility
- Evaporator

Low pressure Switch.



Laboratory Instrument India